


ProteoScan Cancer Lysate Array (PA100002)

Array Layout Scheme

Cancer Lysate Array Layout – PA100002



4.5 mm

Sub-arrays

B1-B3 Breast (B)

C1-C3 Colon (C)

K1-K3 Kidney (K)

L1-L3 Lung (L)

Ly1-Ly3 Lymphoma (Ly)

O1-O3 Ovary (O)

P1-P3 Pancreas (P)


Pr1-Pr3 Prostate (Pr)

Mel1-Mel3 Melanoma (Mel)

Lv1-Lv3 Liver (Lv)

St1-St3 Stomach (St)

Layout - Sub Arrays X1 and X2 – PA100002



BSA-Cy3

BSA-Cy5

BSA-Cy3-Cy5

BSA (ug/ml)

Buffer

Lysate (ug/ml)

Empty

Cy3, Cy5
location
controls

BSA
Background
d controls

Layout - Sub Arrays X3 – PA100002

