

Sure-Race™ Discovery Panels

KIT COMPONENTS

- Four identical, sealed 48-well plates containing two dilutions of 24 human or mouse RACE-ready tissue-specific cDNA's
- Set of adapter primers; ADP 1, outer primer for first-round PCR and ADP 2 for nested, second-round PCR
- Transferrin receptor primers 1 and 2 (10 µM each) positive control primers for first and second-round PCR
- Free sample of OriGene's DNA Quanti-Ladder™, a ready-to-use molecular marker for easy quantitation and size determination (50 µl)

ORDERING INFORMATION & PRICING

SURE-RACE DISCOVERY PANELS	CATALOG NUMBER
Human Sure-RACE™ Panels	HRAA-101
Mouse Sure-RACE™ Panels	MRAB-101

For pricing information call Customer Service at 1-888-267-4436 ext. 1.
To Order call 888-267-4436 (M-F, 9AM-5PM EST) or fax your order to 1-301-340-9254 (24 hours)
For Technical Support call 301-340-3188 ext. 2 or visit the OriGene website at www.origene.com

PRODUCT DELIVERY AND STORAGE

Panels are shipped at room temperature. Long-term storage at -20°C is recommended.



TrueClone™ Collection Blue-Ribbon™ Reagents Rapid-Screen™ Tools

6 Taft Court, Suite 300 Rockville, MD 20850 • Phone: 888-267-4436 Fax: 301-340-9254

www.origene.com